ウエル名	1	2	3	4	5	6	7	8	9
ArgA	v	-	v	v	+	v	-	v	v
GGT	+	-	v	-	v	v	+	v	v
LysA	v	-	+	v	v	-	v	v	v
dGAL	v	-	v	+	v	v	-	v	v
LeuA	v	+	v	+	v	v	v	v	-
ELLM	v	+	v	v	-	v	+	v	v
PheA	v	-	v	+	v	v	+	v	v
ProA	-	+	v	-	v	v	+	v	v
PyrA	v	-	+	-	v	v	v	v	v
TyrA	v	-	v	v	v	v	+	v	v
APPA	-	+	v	-	+	v	v	v	v
dGLU	+	-	v	+	v	v	-	v	v
GLYG	v	-	v	v	v	v	-	+	v
dMNE	+	-	v	v	v	v	-	v	v
dMAL	+	-	v	-	v	v	v	v	v
SAC	+	-	v	v	v	v	-	v	v
NAG	v	-	+	v	v	v	-	v	v
URE	-	-	v	+	v	v	v	v	v
BGALi	+	-	v	-	v	+	v	v	v
ODC	-	+	v	v	-	v	v	v	v
AARA	v	-	v	v	-	+	v	v	v
PVATE	-	-	v	v	v	v	+	v	v
PHC	v	-	v	+	-	v	v	v	v
dMLT	v	v	v	+	-	v	+	v	v
MTE	+	-	v	v	v	v	-	v	v
IGLM	-	v	v	v	-	v	+	v	v
PHOS	+	-	v	+	v	v	-	v	v
dRIB2	v	-	v	+	v	v	-	v	v
OPS	v	-	v	+	v	-	v	v	v
dXYL	v	-	v	+	v	v	-	v	v

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- 2）測定範囲
- 本品で同定される菌種は下記のとおり；
- Actinobacillus ureae*
Campylobacter coli
Campylobacter fetus ssp *fetus*
Campylobacter jejuni ssp *jejuni*
Capnocytophaga spp.
Cardiobacterium hominis
Eikenella corrodens
Gardnerella vaginalis
Haemophilus actinomycetemcomitans
Haemophilus aphrophilus
Haemophilus haemolyticus
Haemophilus influenzae
Haemophilus parahaemolyticus
Haemophilus parainfluenzae
Haemophilus paraphrophilus
Haemophilus segnis
Kingella denitrificans
Kingella kingae
Moraxella catarrhalis
Neisseria cinerea
Neisseria elongata
Neisseria gonorrhoeae
Neisseria lactamica
Neisseria meningitidis

Neisseria sicca
Oligella urethralis
Suttonella indologenes

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- 主要文献
1. Balows A, Hausler WJ, Herrmann KL, Isenberg HD, Shadomy HJ, editors. *Manual of Clinical Microbiology*, 5th Edition. American Society for Microbiology, Washington, D.C. 1991.
2. Balows A, Truper HG, Dworkin M, Harder W, and Schleifer K-H, editors. The Prokaryotes - a Handbook on the Biology of Bacteria: Exophysiology, Isolation, Identification, Applications, 2nd Edition Volume II. Springer-Verlag, New York 1992.
3. Difco Manual Dehydrated Culture Media and Reagents for Microbiology, 10th Edition, 1984.
4. Holt JG, Krieg NR, Sneath PHA, Staley JT, Williams ST, editors. *Bergey's Manual of Determinative Bacteriology*, 9th Edition. William and Wilkins, Baltimore, Maryland. 1994.
5. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC, editors. *Color Atlas and Textbook of Diagnostic Microbiology* 4th Edition. Lippincott, Philadelphia, PA 1992.
6. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC, editors. *Color Atlas and Textbook of Diagnostic Microbiology* 5th Edition. Lippincott, Philadelphia, PA 1997.
7. Krieg NR, Holt JG, editors. *Bergey's Manual of Systematic Bacteriology*, 9th Edition. Williams and Wilkins, Baltimore 1984.
8. Manafi M, Kneifel W, Bascomb S. Fluorogenic and chromogenic substrates used in bacterial diagnostics. Microbiol. 1991; Rev.55:335-348
9. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Yolklen RH, editors. *Manual of Clinical Microbiology*, 7th Edition. American Society for Microbiology, Washington, D.C. 1999.
10. Murray PR, Baron EJ, Jorgensen JH, Pfaller MA and Yolklen RH, editors. *Manual of Clinical Microbiology*, 8th Edition. American Society for Microbiology, Washington, D.C. 2003.

11. National Committee for Clinical Laboratory Standards, M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue* - Approved Guideline, 1997.
12. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, *Biosafety in Microbiological and Biomedical Laboratories*, 1988.
13. Weyant RS, Moss CW, Weaver RE, Hollis DG, Jordon JG, Cook EC, and Daneshvar MI. Identification of Unusual Pathogenic and Gram-Negative Aerobic and Facultatively Anaerobic Bacteria 2nd Edition. Williams & Wilkins, Philadelphia, PA 1996.

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Product covered by one or more of U.S. Patent Numbers D 414,272 ; D 437,797 ; 5,609,828 ; 5,746,980 ; 5,804,437 ; 5,869,005 ; 5,932,177 ; 5,951,952 ; 6,267,929 ; 6,309,890 and 6,340,573 ; and Foreign Counterparts.
Other Patents Pending.

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