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1. AMIES CR. 1967. A modified formula for the preparation of Stuart's medium. Canadian Journal of Public Health, July 1967. Vol. 58, 296-300.
2. HUDSPETH M.K., CITRON D.M., GOLDSTEIN E.J.C. – Evaluation of a novel specimen transport system (Venturi transystem) for anaerobic bacteria. – *Clinical Infectious Disease*, 1997, 25 (suppl 2), S132–S133.
3. ISENBERG HD, WASHINGTON JA II, DOERN GV, AMSTERDAM DN, Specimen Collection and Handling, p. 15-28, Manual of Clinical Microbiology, Fifth Edition, American Society for Microbiology, Washington DC, 1991.
4. ISENBERG HD(Editor in Chief). Clinical Microbiology Procedures handbook. American Society for Microbiology, Washington DC, 1992.
5. PERRY JL. Assessment of swab transport systems for aerobic and anaerobic organism recovery. Journal of Clinical Microbiology, May 1997. Vol. 35, 1269-1271.
6. SUMMANEN P, BARON EJ, CITRON DM, STRONG CA, WEXLER HM, FINEGOLD SM. Wadsworth Anaerobic Bacteriology Manual. Fifth Edition. Star Publishing Company, Belmont CA, 1993.
7. THOMPSON D.S., FRENCH S.A. – Comparison of commercial Amies transport systems with in-house Amies medium for recovery of *Neisseria gonorrhoeae*. – *Journal of Clinical Microbiology*. Vol. 37, n° 9, Sept. 1999, p. 3020–3021.

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