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- 1.Bolmström, A. et al. (1988). A Novel Technique for Direct Quantification of Antimicrobial Susceptibility of Microorganisms. ICAAC, poster 1209.
- 2.Baker, C. N. et al. (1991). Comparison of the Etest to Agar Dilution, Broth Microdilution, and Agar Diffusion Susceptibility Testing Techniques by Using a Special Challenge Set of Bacteria. Journal of Clinical Microbiology.
- 3.Brown, D. F. J. and Brown, L. (1991). Evaluation of the Etest, a novel method of quantifying antimicrobial activity. Journal of Antimicrobial Chemotherapy.
- 4.Jorgensen, J. H. et al. (1994). Detection of penicillin and extended spectrum cephalosporin resistance among *S. pneumoniae* clinical isolates using Etest. Journal of Clinical Microbiology
- 5.Citron D. M. et al. (1991). Evaluation of Etest for susceptibility testing of anaerobic bacteria. Journal of Clinical Microbiology .
- 6.Sanchez M. et al. (1993). Etest, an antimicrobial susceptibility testing method with broad clinical and epidemiological application. The Antimicrobial Newsletter.
- 7.Schulz J. E. et al. (1993). Reliability of the Etest for detection of ampicillin, vancomycin, and high-level aminoglycoside resistance in *Enterococcus* spp. Journal of Clinical Microbiology.
- 8.Baker C. N. et al. (1994). Optimizing testing of methicillin resistant *Staphylococcus* spp. Diagnostic Microbiology and Infectious Disease.
- 9.Tenover F. C. et al. (1996). Evaluation of commercial methods for determining antimicrobial susceptibility of *S. pneumoniae*. Journal of Clinical Microbiology.
- 10.Rosenblatt J. E. et al. (1995). Evaluation of the Etest for susceptibility testing of anaerobic bacteria. Diagnostic Microbiology and Infectious Disease.
- 11.Class II Special Controls Guidance Document: Antimicrobial Susceptibility Test (AST) Systems; Guidance for Industry and FDA. August 2009.
- 12.Lorian, V. Antibiotics in Laboratory Medicine. 5th Ed. 2005. Williams & Wilkins, USA.
- 13.Murray, P.R. et al. Manual of Clinical Microbiology. 9th Ed. 2001. ASM Press.
- 14.CLSI. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically. Approved Standard, M7-A (latest edition).
- 15.CLSI. Methods for Dilution Antimicrobial Susceptibility Tests of Anaerobic Bacteria. Approved Standard, M11-A (latest edition).
- 16.CLSI Performance standards for antimicrobial susceptibility testing. M100-S (latest edition).
- 17.Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard. CLSI M27-A (latest version).
- 18.Pfaller M.A., Messer S.A., Bolmström A., Odds F.C., Rex J.H. (1996). Multisite reproducibility of the Etest method for Antifungal susceptibility testing of yeast isolates. *Journal of Clinical Microbiology (JCM)*. **34**(7):1691-1693.
- 19.Espinel-Ingroff A., Pfaller M.A., Erwin M.E., Jones R.N. (1996). Interlaboratory evaluation of Etest method for testing antifungal susceptibilities of pathogenic yeasts to five antifungal agents by using casitone agar and solidified RPMI 1640 medium with 2% glucose. *JCM*. **34**(4): 848-852.
- 20.Pfaller M.A., Messer S.A., Karlsson Å., Bolmström A. (1998). Evaluation of the Etest method for determining fluconazole susceptibilities of 402 clinical yeast isolates by using three different agar media. *JCM*. **36**(9):2586-2589.
- 21.Pfaller M.A., Messer S.A., Mills K., Bolmström A., Jones R.N. (2000). Evaluation of the Etest method for determining voriconazole susceptibilities of 312 clinical isolates of *Candida* species by using three different agar media. *JCM*. **38**(10): 3715-3717.
- 22.Pfaller M.A., Messer S.A., Mills K., Bolmström A., Jones RN. (2001). Evaluation of the Etest method for determining posaconazole MICs for 314 clinical isolates of *Candida* species. *JCM*. **39**(11): 3952-3954.
- 23.Barry A.L., Pfaller M.A., Rennie R.P., Fuchs P.C., Brown S.D. (2002). Precision and accuracy of fluconazole susceptibility testing by broth microdilution, Etest and disk diffusion methods. *Antimicrobial Agents and Chemotherapy (AAC)*. **46**(6): 1781-1784.
- 24.Maxwell M.J., Messer S.A., Hollis R.J., Boyken L., Tendolkar S., Diekema D.J., Pfaller M.A., and the International Fungal Surveillance Participant Group. (2003). Evaluation of Etest method for determining fluconazole and voriconazole MICs for 279 clinical isolates of *Candida* species infrequently isolated from blood. *JCM*. **41**(3): 1087-1090.
- 25.Pfaller M.A., Diekema D.J., Messer S.A., Boyken L., Hollis R.J. (2003). Activities of fluconazole and voriconazole against 1,586 recent clinical isolates of *Candida* species determined by broth microdilution, disk diffusion and Etest methods: Report from Artemis global antifungal susceptibility program 2001. *JCM*. **41**(4): 1440-1446.
- 26.Pfaller M.A., Diekema D.J., Boyken L., Messer S.A., Tendolkar S., Hollis R.J. (2003). Evaluation of Etest and disk diffusion methods for determining susceptibilities of 235 bloodstream isolates of *C. glabrata* to fluconazole and voriconazole. *JCM*. **41**(5): 1875-1880.
- 27.河喜多龍祥「薬剤感受性検査」 近代出版

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